Paragon

Canada Instructions for Use

Paragon CRT® Contact Lenses, Paragon CRT®100 Contact Lenses and Paragon CRT Dual Axis® Contact Lenses

CAUTIONS:

- Please read carefully and keep this information for future use.
- Contact lenses for Corneal Refractive Therapy should be fitted only by a contact lens fitter trained and certified in the fitting of conventional and sigmoid geometry contact lenses.
- It is essential that you follow your eye care professional's directions and all labeling instructions for proper use of your contact lenses and lens care products.
- Nonsterile. Clean and condition lenses prior to use.

PRODUCT NAMES

Paragon CRT® Contact Lenses

Paragon CRT®100 Contact Lenses

Paragon CRT Dual Axis® Contact Lenses

MATERIAL AND PACKAGING

Paragon CRT® Contact Lenses are manufactured from Paragon HDS® (paflucocon B), Paragon CRT® 100 and Paragon CRT Dual Axis® Contact Lenses are manufactured from Paragon HDS®100 (paflucocon D). These lenses for Corneal Refractive Therapy are to be worn overnight with removal during all or part of each following day. Both materials are thermoset fluorosilicone acrylate copolymer derived primarily from siloxane acrylate, trifluoroethyl methacrylate and methylmethacrylate with a water content of less than 1%. These contact lenses for Corneal Refractive Therapy are available as lathe cut firm contact lenses with blue, red, yellow and green tints. The blue tinted lens contains D&C Green No. 6. The green lens contains D&C Green No. 6 and Perox Yellow No. 9. The red lens contains D&C Red No. 17 and yellow lens contain Perox Yello No. 9. These products may be plasma treated.

Each Paragon CRT® Contact Lens, Paragon CRT®100 and Paragon CRT Dual Axis® Contact Lens is supplied nonsterile in an individual plastic vial. The lens is shipped dry; or, wet shipped in solution. The solution contains poloxamine, hydroxyalkylphosphonate, boric acid, sodium borate, sodium chloride, hydroxypropylmethyl cellulose, Glucam and preserved with chlorhexidine gluconate (0.003%), polyaminopropyl biguanide (0.0005%). The vial labeling is marked with the central base curve radius, dioptric power, overall diameter, Return Zone Depth, Landing Zone Angle, center thickness and serial number.

Never reuse the solution. You may store the lens in the unopened container until ready to dispense, up to the solution expiry date indicated on the vial labeling. If the lens is stored for longer periods of time, it should be cleaned and disinfected with a recommended product. Follow the directions on the selected disinfecting solution regarding prolonged storage.

INTENDED USE

Rigid Gas Permeable Paragon CRT® Contact Lenses (paflucocon B), Paragon CRT®100 and Paragon CRT Dual Axis® Contact Lenses (paflucocon D) for Corneal Refractive Therapy are indicated for use in the reduction of myopic refractive error in nondiseased eyes. The lenses are indicated for overnight wear in a Corneal Refractive Therapy fitting program for the temporary reduction of myopia up to 6.00 diopters in eyes with astigmatism up to 1.75 diopters. The lenses may be disinfected using only a chemical disinfection system.

Note: To maintain the Corneal Refractive Therapy effect of myopia reduction overnight lens wear must be continued on a prescribed schedule. Failure to do so can affect daily activities (e.g., night driving), visual fluctuations and changes in intended correction.

RECOMMENDED WEARING SCHEDULE

Although many practitioners have developed their own initial wearing schedules, the following sequence is recommended as a guideline. Patients are cautioned to limit the wearing schedule recommended by the eye care practitioner regardless of how comfortable the lenses feel. It is ideal for the patient to start with overnight wear the first night. A well fit lens provides for centration on the closed eye. The effects of lid interaction on blinking and gravity may result in lens decentration during open eye wear. Patients should be instructed to place the lens on the eye 15 to 20 minutes before going to sleep.

Patients are cautioned; "when in doubt, take it out". It is important that the new wearer not sleep in a lens that has a significant foreign body sensation. In the event of foreign body sensation, the patient should remove the lens, clean and rewet it and replace the lens. If the sensation continues, the lens should not be worn. The patient should report for follow-up evaluation the morning after the first overnight wear. The visit is best scheduled within a few hours of awakening and the patient should report with the lens in place. This visit provides an excellent opportunity to evaluate lens centration and potential lens adherence.

Upon the absence of clinical signs and complications, the patient may be instructed to continue overnight wear of the lens until the next scheduled follow-up visit.

An alternate initial daytime wear schedule may be offered at the practitioner's discretion.

Day 1 two periods of wear not to exceed 6 hours total

Day 2 6 hours

Day 3-5 8 hours

Day 6 overnight wear with follow-up visit within 24 hours.

The cornea normally changes within five to eight hours of wear. The wearing schedule should be modulated to determine the MINIMUM wear required for myopic reduction. The average wearing time is between 8 and 10 hours. Attempt to maintain wearing time at this level.

With Paragon CRT® Contact Lenses, Paragon CRT® 100 and Paragon CRT Dual Axis® Contact Lenses, the lens used to achieve refractive therapy is usually the lens used to maintain achieved correction. The Retainer Lens wearing time begins with the same wearing time required for the last fitted Paragon CRT® Contact Lenses, Paragon CRT®100 and Paragon CRT Dual Axis®Contact Lenses for Contact Lens Corneal Refractive Therapy. After a period of several days, or when the eye care practitioner is satisfied that the patient has adapted to the first Retainer Lenses, the patient may attempt to skip a night of wear to monitor the duration of visual improvement. This may continue for as long as the patient can see clearly. The schedule of overnight wear must be modulated to maintain visual performance.

HANDLING YOUR LENSES

It is essential that you learn and use good hygienic methods in the care and handling of your new lenses. Cleanliness is the first and most important aspect of proper contact lens care. In particular, your hands should be clean and free of any foreign substance when you handle your lenses. The procedures are:

- Always wash your hands thoroughly with a mild soap, rinse completely, and dry with a lint-free towel before touching your lenses.
- Avoid the use of soaps containing cold cream, lotion, or oily cosmetics before handling your lenses, since these substances may come into contact with the lenses and interfere with successful wearing.
- To avoid damaging your lenses, handle them with your fingertips, and be careful to avoid contact with your fingernails. Do not use tweezers

Inserting The Lens

Work over a table, upon which is placed a clean towel. Do not place lenses on the eye while working over a sink. For the right eye:

- Wet your right index finger with a drop of conditioning solution and place the contact lens front side down on your right index finger.
- Place the second finger of the left hand on the middle of the upper lid and press firmly upward.
- Place the second finger of the right hand on the lower lid and press firmly downward. Stare into a mirror as though looking through the second finger holding the contact lens. You will later learn to do this without a mirror.
- Slowly move the hand to advance the forefinger with the contact lens towards the cornea until the lens touches the cornea and release the lids.
- Release the lid and close the eye for a few seconds.

Repeat procedure for the left eye.

- Always remove the same lens first.
- Work over a table with a clean towel. Do not remove lenses over a sink.
- Place the right index finger at the outer corner of the eye.
- Place the left hand cupped below the eye.
- Open the eyes wide as if to stare.
- Continue to keep the eyes open and pull the lids sideways away from nose.
- Blink quickly and firmly. Remove the second lens by following the same procedure.

OPTIONAL METHOD: If you choose to use a DMV Classic lens remover device, compress the handle of the DMV between your thumb and index finger while looking straight ahead carefully center it squarely on the contact lens. When the Classic has been seated on the lens, release the compression pressure and gently tip the remover down to retract the lens. Always be aware of the position of the lens on the eye before applying the Classic to the lens. Then, squeeze the stem of the remover again to release the suction and slide the lens from the remover between your thumb and forefinger. Remove second lens following the same procedure. Follow the manufacturer's instructions for cleaning your DMV Classic. For all other lens removing devices follow manufacturer's instructions.

Follow the required lens care procedures described under the heading: CARING FOR YOUR LENSES.

Note: If these methods of removing your lens are difficult for you, your eye care practitioner will provide you with an alternate method.

Lens Wearing Precautions

- Regular follow-up visits are necessary to assure the continuing health of your eyes. Adhere to your eye care practitioner's recommended schedule for follow-up visits.
- Never wear your contact lenses beyond the period recommended by your eye care practitioner.
- Always contact your eye care practitioner before using any medicine in your eyes.
- Avoid, if possible, all harmful or irritating vapors and fumes when wearing lenses.
- If aerosol products such as sprays are used while wearing lenses, exercise caution and keep eyes closed until the spray has settled.
- Ask your eye care practitioner about wearing your lenses during sporting activities.
- Inform your employer that you wear contact lenses, especially if your job involves use of eye protection equipment.
- Very rarely, a lens that is on the cornea will be displaced onto the white part of the eye during lens wear. This may also occur during placement and removal of the lenses if the correct techniques are not performed properly. To center a lens follow this procedure.
 - First locate the lens by pulling away the lids.
 - After the lens is found, gently press on the lid over the lens while looking away from the direction of the lens.
 - Next look back towards the lens.
- If the lens sticks (stops moving) or cannot be removed, you should apply 2 drops of the recommended lubricating or rewetting solution directly to the eye and wait until the lens begins to move freely on the eye before removing it. If nonmovement of the lens continues after 30 minutes, you should IMMEDIATELY consult your eye care practitioner or attending hospital emergency room physician.
- WARNING: The risk of ulcerative keratitis has been shown to be greater among wearers of extended wear lenses than among wearers of daily wear lenses. The risk among extended wear lens wearers increases with the number of consecutive days that lenses are worn between removals, beginning with the first overnight use. This risk can be reduced by recognizing the signs of a possible problem (see POSSIBLE PROBLEMS) and by carefully following directions for routine lens care, including cleaning of the lens storage case. Additionally, smoking increases the risk of ulcerative keratitis for contact lens wearers. It is recommended that contact lens wearers see their eye care practitioner twice each year or, if directed, more frequently.

CARING FOR YOUR LENSES

For continued safe and comfortable wearing of your lenses, it is important that you clean and rinse, then disinfect your lenses after each removal using the care regimen recommended by your eye care practitioner. Cleaning and rinsing are necessary to remove mucus, secretions, films, or deposits, which may have accumulated during wearing. The ideal time to clean, rinse and disinfect your lenses is immediately after wearing them. Disinfecting is necessary to destroy harmful germs.

You should adhere to a recommended care regimen. Failure to follow the regimen may result in development of serious ocular complications as discussed in the POSSIBLE PROBLEMS section.

When you first receive your lenses, practice how to put the lenses on and how to remove them while you are in your eye care practitioner's office. At that time you will be provided with a recommended cleaning and disinfection regimen and, instructions and warnings for lens care, handling, cleaning, and disinfection. Your eye care practitioner should instruct you about appropriate and adequate procedures and products for your use.

For safe contact lens wear you should know and always practice your lens care routine.

- Always wash, rinse, and dry hands before handling contact lenses.
- Always use fresh unexpired lens care solutions.
- Use the recommended system of lens care, which is chemical (not heat) and carefully follow instructions in solution labeling. Different solutions cannot always be used together and not all solutions are safe for use with all lenses. **Do not alternate or mix lens care systems unless indicated on solution labeling.**
- Always remove, clean, rinse, enzyme and disinfect your lenses according to the schedule prescribed by your eye care practitioner. The use of an enzyme or any cleaning solution does not substitute for disinfection.
- **To avoid contamination, do not use saliva, tap water or anything other than the recommended solutions for lubricating or rewetting your lenses. Do not put lenses in your mouth.**

For safe contact lens wear you should know and always practice your lens care routine.

- Always wash, rinse, and dry hands before handling contact lenses.
- Always use fresh unexpired lens care solutions
- Use the recommended system of lens care, which is chemical (not heat) and carefully follow instructions in solution labeling. Different solutions cannot always be used together and not all solutions are safe for use with all lenses. Do not alternate or mix lens care systems unless indicated on solution labeling.
- Always remove, clean, rinse, enzyme and disinfect your lenses according to the schedule prescribed by your eye care practitioner. The use of an enzyme or any cleaning solution does not substitute for disinfection.
- To avoid contamination, do not use saliva, tap water or anything other than the recommended solutions for lubricating or rewetting your lenses. Do not put lenses in your mouth.

Always wash and rinse your hands thoroughly before handling your contact lenses.

- Clean one lens first (always start with the same lens first to avoid mix-ups). Place the lens, front side down, in the palm of the hand and apply several drops of cleaning solution. Using the index finger of the other hand, apply slight pressure in a swirling motion for the time recommended by the cleaning solution manufacturer. Do not clean the lens by rubbing it between the thumb and index fingers, as this may cause lens warpage. In addition, your lenses may require periodic enzymatic cleaning to remove protein deposits.
- Rinse the lens thoroughly as recommended by your lens care product manufacturer to remove the cleaning solution, mucus, and film from the lens surface. Place that lens into the correct chamber of the lens storage case. Then repeat the clean and rinse procedure for the second lens.
- After cleaning and rinsing the lenses disinfect them by using the system recommended by your eye care practitioner and/or the lens manufacturer. Follow the instructions provided in the disinfection solution labeling.
- To store lenses, disinfect and leave them in the closed case until ready to wear. If lenses are not to be used immediately following disinfection, you should consult the storage solution package insert or your eye care practitioner for information on storage of your lenses. Always keep your lenses completely immersed in a recommended disinfecting/conditioning solution when the lenses are not being worn. If you discontinue wearing your lenses, but plan to begin wearing them again after a few weeks, ask your eye care practitioner for a recommendation on how to store your lenses.

Note: Paragon CRT® Contact Lenses, Paragon CRT®100 and Paragon CRT Dual Axis® Contact Lenses for Corneal Refractive Therapy cannot be heat (thermally) disinfected.

Contact lens cases can be a source of bacteria growth. After removing your lenses from the lens case, empty and rinse the lens storage case with solution(s) recommended by the lens case manufacturer; then allow the lens case to air dry. When the case is used again, refill it with fresh disinfecting solution. Lens cases should be replaced at regular intervals as recommended by the lens case manufacturer or your eye care practitioner.

WHEN NOT TO WEAR LENSES

DO NOT USE Paragon CRT®, Paragon CRT®100 or Paragon CRT Dual Axis® Contact Lenses for Corneal

Refractive Therapy when any of the following conditions exist:

- Acute and subacute inflammations or infection of the anterior segment of the eye.
- Any eye disease, injury, or abnormality that affects the cornea, conjunctiva or eyelids.
- Severe insufficiency of tears (dry eyes).
- Corneal hypoesthesia (reduced corneal sensitivity).
- Any systemic disease which may affect the eye or be exacerbated by wearing contact lenses.
- Allergic reactions of ocular surfaces or adnexa which may be induced or exaggerated by wearing contact lenses or use of contact lens solutions.
- Allergy to any ingredient, such as mercury or thimerosal, in a solution which is to be used to care for contact lenses.
- Any active corneal infection (bacterial, fungal or viral).
- If eyes become red or irritated.

POSSIBLE PROBLEMS

Although contact lenses provide many benefits to the wearer, it is possible that problems might occur and may be first noticed as one or more of the following conditions:

- Eyes stinging, burning, itching (irritation), or other eye pain
- Comfort is less than when lens was first placed on eye
- Feeling of something in the eye such as a foreign body or scratched area
- Excessive watering (tearing) of the eyes
- Unusual eye secretions
- Redness of the eyes
- Reduced sharpness of vision (poor visual acuity)
- Blurred vision, rainbows, or halos around objects
- Sensitivity to light (photophobia)
- Dry eyes

If the patient notices any of these conditions, the patient should IMMEDIATELY REMOVE THE LENSES. The patient should follow these instructions.

- If the discomfort or problem stops, then look closely at the lens. If the lens is in any way damaged, DO NOT put the lens back on your eye. Place the lens in the storage case and contact your eye care practitioner.
- If the lens has dirt, an eyelash, or other foreign objects on it, or the problem stops and the lens appears undamaged, you should thoroughly clean, rinse and disinfect the lens; then reinsert it.
- If the problem continues, you should IMMEDIATELY remove the contact lenses and consult your eye care practitioner.

When any of the above problems occurs, a serious condition such as infection, corneal ulcer, neovascularization, iritis, persistent stromal edema or GPC (giant papillary conjunctivitis) may be present. Keep the lens off your eye and seek immediate professional identification of the problem and prompt treatment to avoid serious eye damage, including corneal scarring, opacification, blindness or loss of eye.

Manufacturer:
Paragon Vision Sciences, Inc
2120 W. Guadalupe Rd.
Gilbert, AZ 85233 USA
1800-528-8279

Lens Specifications Abbreviations			
Power	PWR	Return Zone Depth	RZD
Base Curve	BC	Landing Zone Angle	LZA
Center Thickness	CT	Overall Diameter	OAD

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